



*Providing Clinical Confidence for Line Placement.*

**AccuSITE™**   
Guiding Patient Care.

1

## attach bracket

Attach custom reusable non-sterile bracket to transducer by aligning locating features.



2

## cover up

CIVCO makes it easier to reduce infection rates and cross-contamination with a telescopically-folded CIV-Flex™ cover. Expand your sterile field to 58" by covering the transducer and cord with the sterile latex-free cover.


safe. efficient. g


Five steps for clinical c

## attach guide

Using proper sterile technique, hook needle guide onto bracket and snap into place with a single click.

4






## depth selection

CIVCO has easy-to-identify color coded needle guides available in 18, 20 and 21 gauge. Following assessment of the vessel, select appropriate depth (0.5, 1.0, 1.5, 2.0, 2.5, 3.0 and 3.5cm) to access targeted vessel.

3

# guiding patient care.

confidence in central line placement.



5

## quick-release

Depress tab on needle guide toward bracket to activate quick-release. Remove needle from transducer assembly.

# AccuSITE™ Transverse Needle Guide



Ref. #	Product Description	Qty.
--------	---------------------	------

## Needle Guide Replacement Kits

610-1017	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-18GA only (Pink)	24
610-1018	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-20GA only (Yellow)	24
610-1019	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-21/22GA only (Green)	24

## Transducer Covers

610-1000	Sterile 10.2 x 147cm (4" x 58") CIV-Flex cover	24
610-989	Sterile 15.2 x 147cm (6" x 58") telescopically-folded polyethylene cover	24

Now available for leading equipment manufacturers including: GE Healthcare, TeleFlex, Terason, Ultrasonix and ZONARE.

To Order Call: 800.445.6741 • 319.656.4447

**AccuSITE™**   
Guiding Patient Care.

# Providing Clinical Confidence for Central Line Placement.

## Safe. Efficient. Guiding Patient Care.

Introducing AccuSITE™, a needle guide designed specifically for guidance in central line placement procedures.

Designed with minimizing hospital risks and increasing patient well being in mind, the AccuSITE needle guide increases first-time success rates during catheter placement procedures.

## Real-time Guidance

Use of real-time and portable ultrasound imaging is increasing the safety and efficiency of infusion therapy. The AccuSITE needle guide coupled with the use of ultrasound has increased first attempt insertion success rates resulting in clinical confidence for central line placement. A recent study shows that with ultrasound guidance, the rate of successful bedside insertions increased to more than 90% compared to a 73% success rate without ultrasound guidance.<sup>1,2</sup>

## Accurate Placement

Consisting of a two-part system of a reusable bracket and a disposable snap-on needle guide, the AccuSITE provides a safe, reliable and efficient tool for vessel cannulation. The AccuSITE holds the needle as it is advanced, providing a smooth insertion, while allowing you to see what you puncture as well as reduce complications by avoiding puncturing adjacent structures.

## Needle Guide Features:


- Single-use design reduces risks of cross-contamination
- Snaps onto bracket with single click- no additional assembly required
- Large colored tab for simple quick-release function
- Large funnel for instrument insertion point

## Transducer Cover Features:

- Improves infection control rates by preventing cross-contamination at minimal cost per procedure
- Longest length cover on market for complete sterile coverage of ultrasound transducer and cord
- Telescopically-folded for easy single-person application to equipment
- Latex-free kit includes sterile gel and bands

## INCREASED SUCCESS RATE<sup>2</sup>

Percentage of Patients



1- "Improved care and reduced costs for patients requiring peripherally inserted central catheters: the role of the bedside ultrasound and a dedicated team" JPEN: Journal of Enteral Nutrition, September/October 2005, p.1.

2- Moureau, Nancy, CRNI, BSN. "Using Ultrasound to GUIDE PICC Insertion." Nursing, December 2003, p. 1.

*Providing Clinical Confidence for Line Placement.*

# AccuSITE™



Guiding Patient Care.

***Placing an Order***

Contact a CIVCO customer service specialist to place an order or request a product sample.  
For international distribution, visit CIVCO Worldwide at [WWW.CIVCO.COM](http://WWW.CIVCO.COM)

Toll-free within the U.S. and Canada: 800.445.6741

Toll-free Facsimile: 877.329.2482

Direct Dial: 319.656.4447

Facsimile: 319.656.4451

E-mail: [info@civco.com](mailto:info@civco.com)

[WWW.CIVCO.COM](http://WWW.CIVCO.COM)

